


Phone Page a.k.a. Conversation Page for Vanguard

© 2003 Musselwhite & Wagner


Why? The phone page is set up for quick conversation while talking on the telephone. Many of the same features (greetings / closings / topic setters) found on the phone page are helpful for conversations, especially for individuals who are switch users. This was set up for a 14-year-old who has accurate but fairly slow access . . . and lots to say!

How? Here are the steps we used to make conversations easier, while still allowing the student to access Minspeak, and allow for quick responses.

Step 1 - Go to Phone Page: First, we started on the basic Minspeak page, then moved to the basic Phone Page, (on R.)



Step 2: Rearrange symbols (using Swap Page Keys from the Toolbox)
Create a row of seven blank symbols, to be used to enter messages for an *interview* (Hey, may I ask you a few questions / first, how many helpings of fruit do you eat each day), *story* (ex: fell in the pool / dog got lost / went to ball game), *directions*, etc. Add **interjection button**. **Be sure to add the tool <Set page temporary> and link it to the Interjection Page (Step 6).**



Step 3: Add #s 1 - 7 at the on the blank row. This will help users move from step to step in the story, but will not require users or helpers to add new symbols to represent each story. If desired, add symbols for messages such as "Guess what" or "I'm not telling."

